

Anti-ADAMTS1 Mouse Monoclonal Antibody [Clone ID: OTI4H2]

Catalog Number: M03009


Overview

Product Name	Anti-ADAMTS1 Mouse Monoclonal Antibody [Clone ID: OTI4H2]
Reactive Species	Human
Description	Boster Bio ADAMTS1 mouse monoclonal antibody, clone OTI4H2 (formerly 4H2). Catalog# M03009. Tested in IHC, WB. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal OTI4H2
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9UHI8


Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 412-644 of human ADAMTS1 (NP_008919) produced in E.coli.
Isotype	IgG2a
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB: 1:2000 IHC: 1:150

Anti-ADAMTS1 Mouse Monoclonal Antibody [Clone ID: OTI4H2] (M03009) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ADAMTS1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ADAMTS1.



Immunohistochemical staining of paraffin-embedded Carcinoma of Human bladder tissue using anti-ADAMTS1 mouse monoclonal antibody. (M03009) Dilution: 1:150

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ADAMTS1 Mouse Monoclonal Antibody [Clone ID: OTI4H2]

For Research Use Only. Not for use in diagnostic procedures.