

Anti-ERK 3 MAPK6 Monoclonal Antibody

Catalog Number: M03011


Overview

Product Name	Anti-ERK 3 MAPK6 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-ERK 3 MAPK6 Monoclonal Antibody catalog # M03011. Tested in ELISA, Flow Cytometry, WB applications. This antibody reacts with Human.
Application	ELISA, Flow Cytometry, WB
Clonality	Monoclonal 400000000
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q16659


Technical Details

Immunogen	Purified recombinant fragment of human ERK 3 expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 FC 1:200-1:400 ELISA 1:10000

Anti-ERK 3 MAPK6 Monoclonal Antibody (M03011) Images



Western Blot analysis using ERK 3 Monoclonal Antibody against HEK293 (1) and ERK3-hlgGfc transfected HEK293 (2) cell lysate.



Flow cytometric analysis of HeLa cells using ERK 3 Monoclonal Antibody (green) and negative control (purple).

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ERK 3 MAPK6 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.