

## Anti-Calcineurin A (PPP3CA) Mouse Monoclonal Antibody [Clone ID: OTI5C6]

Catalog Number: M03026-1


### Overview

Product Name	Anti-Calcineurin A (PPP3CA) Mouse Monoclonal Antibody [Clone ID: OTI5C6]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio PPP3CA mouse monoclonal antibody, clone OTI5C6. Catalog# M03026-1. Tested in IHC, WB. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Monoclonal OTI5C6
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q08209


### Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 2-105 of human PPP3CA (NP_000935) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500 IHC 1:2000


## Anti-Calcineurin A (PPP3CA) Mouse Monoclonal Antibody [Clone ID: OTI5C6] (M03026-1) Images



Western blot analysis of extracts (35ug) from cell lines by using anti-PPP3CA monoclonal antibody (1:500).



Immunohistochemical staining of paraffin-embedded Human adult brain tissue within the normal limits using anti-PPP3CA mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)



Immunohistochemical staining of paraffin-embedded Human embryonic cerebellum within the normal limits using anti-PPP3CA mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Calcineurin A (PPP3CA) Mouse Monoclonal Antibody [Clone ID: OTI5C6]

For Research Use Only. Not for use in diagnostic procedures.