

## Anti-Collagen XVII Rabbit Monoclonal Antibody

Catalog Number: M03031


### Overview

Product Name	Anti-Collagen XVII Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Collagen XVII Rabbit Monoclonal Antibody catalog # M03031. Tested in WB, IHC, ICC/IF, IP applications. This antibody reacts with Human, Mouse, Rat.
Application	IP, IF, IHC, ICC, WB
Clonality	Monoclonal AA-3
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UMD9


### Technical Details

Immunogen	A synthesized peptide derived from human Collagen XVII
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 ICC/IF 1:50-200 IP 1:30


## Anti-Collagen XVII Rabbit Monoclonal Antibody (M03031) Images




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Immunohistochemical analysis of paraffin-embedded mouse skin, using Collagen XVII Antibody.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Immunofluorescent analysis of A431 cells, using Collagen XVII Antibody .

Western blot analysis of Collagen XVII expression in Human skin lysate.



## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Collagen XVII Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.