

## Anti-UBE2S Mouse Monoclonal Antibody [Clone ID: OTI4H3]

Catalog Number: M03045


### Overview

Product Name	Anti-UBE2S Mouse Monoclonal Antibody [Clone ID: OTI4H3]
Reactive Species	Human, Monkey, Mouse, Rat
Description	Boster Bio UBE2S mouse monoclonal antibody, clone OTI4H3 (formerly 4H3). Catalog# M03045. Tested in FC, IF, WB. This antibody reacts with Human, Monkey, Mouse, Rat.
Application	Flow Cytometry, IF, WB
Clonality	Monoclonal OTI4H3
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q16763


### Technical Details

Immunogen	Full length human recombinant protein of human UBE2S (NP_055316) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:200~2000 IF 1:100 Flow Cytometry 1:100


## Anti-UBE2S Mouse Monoclonal Antibody [Clone ID: OTI4H3] (M03045) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY UBE2S (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UBE2S.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-UBE2S monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).



HEK293T cells transfected with either UBE2S (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-UBE2S antibody (M03045)



Anti-UBE2S mouse monoclonal antibody (M03045) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY UBE2S.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-UBE2S Mouse Monoclonal Antibody [Clone ID: OTI4H3]

For Research Use Only. Not for use in diagnostic procedures.