

## Anti-Calbindin-D Calb1 Antibody (Monoclonal, 31C05)

Catalog Number: M03047-3

### About Calb1

The human calbindin D28k (CALB1) is a calcium-binding protein belonging to the troponin C superfamily. It was originally described as a 27-kD protein induced by vitamin D in the duodenum of the chick. Parmentier et al. (1991) mapped the CALB1 gene to 8q21.3-q22.1 by in situ hybridization.

### Overview

Product Name	Anti-Calbindin-D Calb1 Antibody (Monoclonal, 31C05)
Reactive Species	Bovine, Human, Mouse, Rat
Description	Boster Bio Anti-Calbindin-D Calb1 Antibody (Monoclonal, 31C05) catalog # M03047-3. Tested in IF, IHC, WB applications. This antibody reacts with Bovine, Human, Mouse, Rat.
Application	IF, IHC, WB
Clonality	Monoclonal 31C05
Formulation	Mouse IgG in stabilizing components, 1.2% sodium acetate and 0.01mg NaN <sub>3</sub> . This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P07171

### Technical Details

Immunogen	Bovine kidney calbindin-D
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Mouse IgG (EK1001) for Western blot, and HRP Conjugated anti-Mouse IgG Super Vision Assay Kit (SV0001-1) for IHC(P).
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Mouse IgG1
Form	Lyophilized
Concentration	Adding 1 ml of PBS buffer will yield a concentration of 100 ug/ml.

Purification	Ascites
Suggested Dilutions	Immunohistochemistry (Paraffin-embedded Section), 1-2ug/ml, Human, bovine, mouse, rat Western blot, 0.5-1ug/ml, Human, bovine, mouse, rat Immunofluorescence, 2ug/ml, Mouse

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Calbindin-D Calb1 Antibody (Monoclonal, 31C05)

For Research Use Only. Not for use in diagnostic procedures.