

Anti-TRAF5 Mouse Monoclonal Antibody [Clone ID: OT17A3]

Catalog Number: M03076


Overview

Product Name	Anti-TRAF5 Mouse Monoclonal Antibody [Clone ID: OT17A3]
Reactive Species	Human
Description	Boster Bio TRAF5 mouse monoclonal antibody, clone OT17A3. Catalog# M03076. Tested in WB. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal OT17A3
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	O00463


Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 334-557 of human TRAF5 (NP_004610) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500

Anti-TRAF5 Mouse Monoclonal Antibody [Clone ID: OTI7A3] (M03076) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TRAF5 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TRAF5 (1:500).



Western blot analysis of extracts (35ug) from 6 different cell lines by using anti-TRAF5 monoclonal antibody (1:500).

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRAF5 Mouse Monoclonal Antibody [Clone ID: OTI7A3]

For Research Use Only. Not for use in diagnostic procedures.