

Anti-Munc 13 4 (UNC13D) Mouse Monoclonal Antibody [Clone ID: OT1F4]

Catalog Number: M03088


Overview

Product Name	Anti-Munc 13 4 (UNC13D) Mouse Monoclonal Antibody [Clone ID: OT1F4]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio UNC13D mouse monoclonal antibody, clone OT1F4 (formerly 1F4). Catalog# M03088. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OT1F4
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q70J99

Technical Details

Immunogen	Full length human recombinant protein of human UNC13D (NP_954712) produced in HEK293T cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:200~1000

Anti-Munc 13 4 (UNC13D) Mouse Monoclonal Antibody [Clone ID: OTI1F4] (M03088) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY UNC13D (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UNC13D.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Munc 13 4 (UNC13D) Mouse Monoclonal Antibody [Clone ID: OTI1F4]

For Research Use Only. Not for use in diagnostic procedures.