

## Anti-DSG3 Reference Antibody (Forerunner patent anti-DSG3)

Catalog Number: M03093-2


### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-DSG3 Reference Antibody (Forerunner patent anti-DSG3)   |
| Reactive Species     | Human  |
| Description          | Boster Bio Anti-DSG3 Reference Antibody (Forerunner patent anti-DSG3) (Catalog # M03093-2). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.774EU/ug,determined by LAL method. Expression system: CHO Cell |
| Application          | ELISA, Flow Cytometry, Functional Assay, Kinetics  |
| Clonality            | Monoclonal   |
| Formulation          | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.  |
| Storage Instructions | Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.  |
| Uniprot ID           | P32926   |


### Technical Details

|         |        |
|---------|--------|
| Isotype | IgG1   |
| Form    | Liquid |

## Anti-DSG3 Reference Antibody (Forerunner patent anti-DSG3) (M03093-2) Images



Anti-DSG3 Reference Antibody (Forerunner patent anti-DSG3) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-DSG3 Reference Antibody (Forerunner patent anti-DSG3) is more than 99.63%

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DSG3 Reference Antibody (Forerunner patent anti-DSG3)

For Research Use Only. Not for use in diagnostic procedures.