

## Anti-TNFRSF18 Mouse Monoclonal Antibody [Clone ID: OTI9G8]

Catalog Number: M03125-2


### Overview

Product Name	Anti-TNFRSF18 Mouse Monoclonal Antibody [Clone ID: OTI9G8]
Reactive Species	Human
Description	Boster Bio TNFRSF18 mouse monoclonal antibody, clone OTI9G8 (formerly 9G8). Catalog# M03125-2. Tested in IHC, WB. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal OTI9G8
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9Y5U5


### Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 26-162 of human TNFRSF18 (NP_004186) produced in E.coli.
Isotype	IgG2b
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IHC 1:150


## Anti-TNFRSF18 Mouse Monoclonal Antibody [Clone ID: OTI9G8] (M03125-2) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TNFRSF18 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TNFRSF18.



Immunohistochemical staining of paraffin-embedded Human lymphoma tissue using anti-TNFRSF18 mouse monoclonal antibody. (M03125-2; heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)



Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human colon tissue using anti-TNFRSF18 mouse monoclonal antibody. (M03125-2; heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNFRSF18 Mouse Monoclonal Antibody [Clone ID: OTI9G8]

For Research Use Only. Not for use in diagnostic procedures.