

## Anti-TNFRSF18 / GITR / CD357 Reference Antibody (ragifilimab)

Catalog Number: M03125-4


### Overview

Product Name	Anti-TNFRSF18 / GITR / CD357 Reference Antibody (ragifilimab)
Reactive Species	Cynomolgus monkey, Human
Description	Boster Bio Anti-TNFRSF18 / GITR / CD357 Reference Antibody (ragifilimab) (Catalog # M03125-4). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human. Endotoxin: < 0.775EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q9Y5U5


### Technical Details

Isotype	IgG1
Form	Liquid

## Anti-TNFRSF18 / GITR / CD357 Reference Antibody (ragifilimab) (M03125-4) Images



Anti-TNFRSF18/GITR/CD357 Reference Antibody (ragifilimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-TNFRSF18/GITR/CD357 Reference Antibody (ragifilimab) is more than 100%

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNFRSF18 / GITR / CD357 Reference Antibody (ragifilimab)

For Research Use Only. Not for use in diagnostic procedures.