

Anti-TRIM33 (C-terminus) monoclonal Antibody

Catalog Number: M03133


Overview

Product Name	Anti-TRIM33 (C-terminus) monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-TRIM33 (C-terminus) monoclonal Antibody catalog # M03133. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 3A6
Formulation	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9UPN9

Technical Details

Immunogen	Synthetic peptide, corresponding to the C-terminus of Human TRIM33.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000

Anti-TRIM33 (C-terminus) monoclonal Antibody (M03133) Images



Western blot detection of TRIM33 (C-terminus) antibody in MCF7&HeLaNE cells lysates using TRIM33 (C-terminus) antibody (1:1000 diluted) .

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRIM33 (C-terminus) monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.