

Anti-Clathrin Heavy Chain 1 Rabbit Monoclonal Antibody

Catalog Number: M03134-4

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Clathrin Heavy Chain 1 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of Clathrin Heavy Chain 1 protein. Validated for research with WB,IHC,IF,ELISA applications. Anti Clathrin Heavy Chain 1 is reactive for Human, Mouse, Rat samples
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	Q00610


Technical Details

Isotype	IgG, Kappa
Purification	Protein A
Suggested Dilutions	IHC 1:50-1:200 WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200


Anti-Clathrin Heavy Chain 1 Rabbit Monoclonal Antibody (M03134-4) Images




Mouse colon was stained with anti-Clathrin Heavy Chain 1 rabbit antibody




Human colon carcinoma was stained with anti-Clathrin Heavy Chain 1 rabbit antibody




Human kidney was stained with anti-Clathrin Heavy Chain 1 rabbit antibody



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-Clathrin Heavy Chain antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: U-251 MG Predicted band size: 192kDa Observed band size: 180kDa



Various whole cell lysates were separated by 4-8% SDS-PAGE, and the membrane was blotted with anti-Clathrin Heavy Chain 1 antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: HeLa Lane 2: Mouse brain Predicted band size: 192kDa Observed band size: 180kDa



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-Clathrin Heavy Chain 1 antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1:Rat brain Lane 2: U-251 MG Predicted band size: 192kDa Observed band size: 180kDa

[Anti-Clathrin Heavy Chain 1 Rabbit Monoclonal Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.