

Anti-GSK alpha GSK3A Monoclonal Antibody

Catalog Number: M03152

Overview

Product Name	Anti-GSK alpha GSK3A Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-GSK alpha GSK3A Monoclonal Antibody catalog # M03152. Tested in ELISA, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 6G2F1 / 6G2F2
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P49840


Technical Details

Immunogen	Purified recombinant fragment of GSK3alpha expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

Anti-GSK alpha GSK3A Monoclonal Antibody (M03152) Images



Western Blot analysis using GSK3alpha Monoclonal Antibody against truncated GSK3 alpha recombinant protein (1) and A549 cell lysate (2).



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma (left) and skin carcinoma (right), showing cytoplasmic localization with DAB staining using GSK3alpha Monoclonal Antibody.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GSK alpha GSK3A Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.