

Anti-Cardiac Troponin C Rabbit Monoclonal Antibody

Catalog Number: M03153-1

About TNNC1

Receptor for bradykinin. It is associated with G proteins that activate a phosphatidylinositol-calcium second messenger system.


Overview

Product Name	Anti-Cardiac Troponin C Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Cardiac Troponin C Rabbit Monoclonal Antibody catalog # M03153-1. Tested in WB, IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Monoclonal 20T71
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P63316


Technical Details

Immunogen	A synthesized peptide derived from human Cardiac Troponin C
Isotype	IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used: WB 1:500-1:2000 IHC 1:50-1:200

Anti-Cardiac Troponin C Rabbit Monoclonal Antibody (M03153-1) Images



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of Cardiac Troponin C expression in Human fetal heart lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cardiac Troponin C Rabbit Monoclonal Antibody