

Anti-IFNGR2 Mouse Monoclonal Antibody [Clone ID: OTI1C2]

Catalog Number: M03158


Overview

Product Name	Anti-IFNGR2 Mouse Monoclonal Antibody [Clone ID: OTI1C2]
Reactive Species	Human
Description	Boster Bio IFNGR2 mouse monoclonal antibody, clone OTI1C2 (formerly 1C2). Catalog# M03158. Tested in IHC, WB. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal OTI1C2
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P38484


Technical Details

Immunogen	Full length human recombinant protein of human IFNGR2 (NP_005525) produced in HEK293T cell.
Isotype	IgG2b
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500 IHC 1:150

Anti-IFNGR2 Mouse Monoclonal Antibody [Clone ID: OTI1C2] (M03158) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY IFNGR2 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-IFNGR2.



Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human colon tissue using anti-IFNGR2 mouse monoclonal antibody (M03158) at 1:150 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IFNGR2 Mouse Monoclonal Antibody [Clone ID: OTI1C2]

For Research Use Only. Not for use in diagnostic procedures.