

Anti-DR3/TNFRSF25 Antibody (Monoclonal, 31T47)

Catalog Number: M03227

About TNFRSF25

TNFRSF25 (Tumor Necrosis Factor Receptor Superfamily Member 25), also known as LARD, APO3, DR3 or TNFR25, is a protein that in humans is encoded by the TNFRSF25 gene. Members of the mammalian tumor necrosis factor receptor (TNFR) family are cell-surface proteins that interact with a corresponding TNF-related ligand family. By fluorescence in situ hybridization, Marsters et al. (1996) mapped the Apo3 gene to 1p36.3. Marsters et al. (1996) showed that ectopic expression of Apo3 in mammalian cells triggered apoptosis and activated the transcription factor NF-kappa-B. They suggested that, like TNFR1, Apo3 may regulate distinct signaling pathways in different cellular contexts.

Overview

Product Name	Anti-DR3/TNFRSF25 Antibody (Monoclonal, 31T47)
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DR3/TNFRSF25 Antibody (Monoclonal, 31T47) catalog # M03227. Tested in WB, ICC/IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, ICC, WB
Clonality	Monoclonal 31T47
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q93038

Technical Details

Immunogen	Recombinant protein within human DR3 / LARD aa 3-197.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunocytochemistry/Immunofluorescence, 1:50-200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DR3/TNFRSF25 Antibody (Monoclonal, 31T47)

For Research Use Only. Not for use in diagnostic procedures.