

Anti-Nucleostemin (GNL3) Mouse Monoclonal Antibody [Clone ID: OTI6F4]

Catalog Number: M03283


Overview

| | |
|----------------------|--|
| Product Name | Anti-Nucleostemin (GNL3) Mouse Monoclonal Antibody [Clone ID: OTI6F4] |
| Reactive Species | Human |
| Description | Boster Bio GNL3 mouse monoclonal antibody, clone OTI6F4. Catalog# M03283. Tested in WB. This antibody reacts with Human. |
| Conjugate | Unconjugated |
| Application | WB |
| Clonality | Monoclonal OTI6F4 |
| Formulation | PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Mouse |
| Uniprot ID | Q9BVP2 |


Technical Details

| | |
|---------------------|---|
| Immunogen | Full length human recombinant protein of human GNL3 (NP_055181) produced in E.coli. |
| Isotype | IgG1 |
| Concentration | 1 mg/ml |
| Purification | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Suggested Dilutions | WB 1:500 |

Anti-Nucleostemin (GNL3) Mouse Monoclonal Antibody [Clone ID: OTI6F4] (M03283) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY GNL3 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-GNL3 (1:500).



Western blot analysis of extracts (35ug) from 4 cell lines lysates by using anti-GNL3 monoclonal antibody (1:500).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Nucleostemin (GNL3) Mouse Monoclonal Antibody [Clone ID: OTI6F4]

For Research Use Only. Not for use in diagnostic procedures.