

## Anti-p18 INK4c (CDKN2C) Mouse Monoclonal Antibody [Clone ID: OTI6D3]

Catalog Number: M03299


### Overview

Product Name	Anti-p18 INK4c (CDKN2C) Mouse Monoclonal Antibody [Clone ID: OTI6D3]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio CDKN2C (p18 INK4c) mouse monoclonal antibody, clone OTI6D3 (formerly 6D3). Catalog# M03299. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI6D3
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P42773


### Technical Details

Immunogen	Full length human recombinant protein of human CDKN2C (NP_001253) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:  For protocols, please visit <a href="https://www.bosterbio.com/protocol-and-troubleshooting/">https://www.bosterbio.com/protocol-and-troubleshooting/</a>

## Anti-p18 INK4c (CDKN2C) Mouse Monoclonal Antibody [Clone ID: OTI6D3] (M03299) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CDKN2C (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CDKN2C.



Equivalent amounts of cell lysates (10 ug per lane) of wild-type HeLa cells (WT) and CDKN2C-Knockout HeLa cells (KO) were separated by SDS-PAGE and immunoblotted with anti-CDKN2C monoclonal antibody M03299

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p18 INK4c (CDKN2C) Mouse Monoclonal Antibody [Clone ID: OTI6D3]