

Anti-ASNS/Asparagine Synthetase Rabbit Monoclonal Antibody

Catalog Number: M03302

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-ASNS/Asparagine Synthetase Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ASNS/Asparagine Synthetase Rabbit Monoclonal Antibody catalog # M03302. Tested in WB, IHC, ICC/IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, ICC, WB
Clonality	Monoclonal FGH-1
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P08243

Technical Details


Immunogen	A synthesized peptide derived from human ASNS
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:1000-1:2000

	IHC 1:50-1:200 ICC/IF 1:50-1:200
--	-------------------------------------

Anti-ASNS/Asparagine Synthetase Rabbit Monoclonal Antibody (M03302) Images



Western blot analysis of ASNS expression in K562 cell lysate.



Immunohistochemical analysis of paraffin-embedded human bladder, using ASNS Antibody.

Anti-ASNS/Asparagine Synthetase Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.