

Anti-GCGR Reference Antibody (crotedumab)

Catalog Number: M03342

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-GCGR Reference Antibody (crotedumab)
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-GCGR Reference Antibody (crotedumab) (Catalog # M03342). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Rat, Human, Mouse. Endotoxin: < 0.793EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P47871


Technical Details

Isotype	IgG4SP
Form	Liquid


Anti-GCGR Reference Antibody (crotedumab) (M03342) Images



Anti-GCGR Reference Antibody (crotedumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-GCGR Reference Antibody (crotedumab) is more than 96.61%



Human GPCR His EGFP HEK293 cells were stained with Anti-GCGR Reference Antibody (crotedumab) and negative control protein respectively

Anti-GCGR Reference Antibody (crotedumab)

For Research Use Only. Not for use in diagnostic procedures.