

Anti-PI3 Kinase p85 beta PIK3R2 Monoclonal Antibody

Catalog Number: M03363-1

Overview

| | |
|----------------------|---|
| Product Name | Anti-PI3 Kinase p85 beta PIK3R2 Monoclonal Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-PI3 Kinase p85 beta PIK3R2 Monoclonal Antibody catalog # M03363-1. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse, Rat. |
| Application | WB |
| Clonality | Monoclonal 8D9-D5-F8 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Mouse |
| Uniprot ID | O00459 |

Technical Details

| | |
|---------------------|---|
| Immunogen | Purified recombinant human PI3 Kinase p85 beta protein fragments expressed in E.coli. |
| Cross Reactivity | No cross reactivity with other proteins. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Suggested Dilutions | WB 1:1000 |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PI3 Kinase p85 beta PIK3R2 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.