

Anti-NOB1 Mouse Monoclonal Antibody [Clone ID: OTI1C12]

Catalog Number: M03504


Overview

Product Name	Anti-NOB1 Mouse Monoclonal Antibody [Clone ID: OTI1C12]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio NOB1 mouse monoclonal antibody, clone OTI1C12 (formerly 1C12). Catalog# M03504. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI1C12
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9ULX3

Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 2-234 of human NOB1 (NP_054781) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000

Anti-NOB1 Mouse Monoclonal Antibody [Clone ID: OTI1C12] (M03504) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NOB1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NOB1 (1:2000).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NOB1 Mouse Monoclonal Antibody [Clone ID: OTI1C12]

For Research Use Only. Not for use in diagnostic procedures.