

Anti-SOST / Sclerostin Reference Antibody (blosozumab)

Catalog Number: M03551


Overview

| | |
|----------------------|--|
| Product Name | Anti-SOST / Sclerostin Reference Antibody (blosozumab) |
| Reactive Species | Cynomolgus monkey, Human, Rat |
| Description | Boster Bio Anti-SOST / Sclerostin Reference Antibody (blosozumab) (Catalog # M03551). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Rat, Cynomolgus monkey, Human. Endotoxin: < 0.799EU/ug,determined by LAL method. Expression system: CHO Cell |
| Application | ELISA, Flow Cytometry, Functional Assay, Kinetics |
| Clonality | Monoclonal |
| Formulation | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column. |
| Storage Instructions | Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles. |
| Uniprot ID | Q9BQB4 |


Technical Details

| | |
|---------|--------|
| Isotype | IgG4SP |
| Form | Liquid |

Anti-SOST / Sclerostin Reference Antibody (blosozumab) (M03551) Images



Anti-SOST/Sclerostin Reference Antibody (blosozumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-SOST/Sclerostin Reference Antibody (blosozumab) is more than 98.7%

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SOST / Sclerostin Reference Antibody (blosozumab)

For Research Use Only. Not for use in diagnostic procedures.