

## Anti-ANP32A Monoclonal Antibody

Catalog Number: M03625


### Overview

Product Name	Anti-ANP32A Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-ANP32A Monoclonal Antibody catalog # M03625. Tested in ELISA, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 8D2
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P39687


### Technical Details

Immunogen	Purified recombinant fragment of human ANP32A expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

## Anti-ANP32A Monoclonal Antibody (M03625) Images



Western blotting validation for Anti-ANP32A Monoclonal Antibody M03625 Western Blot (WB) analysis using ANP32A Monoclonal Antibody against HeLa (1) cell lysate. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration



Immunohistochemistry validation of ANP32A using Anti-ANP32A Monoclonal Antibody (M03625). Immunohistochemistry (IHC) analysis of paraffin-embedded Human Kidney tissue with DAB staining using ANP32A Monoclonal Antibody. For more protocol information of IHC

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-ANP32A Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.