

Anti-PAR4 PAWR Monoclonal Antibody

Catalog Number: M03637


Overview

Product Name	Anti-PAR4 PAWR Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-PAR4 PAWR Monoclonal Antibody catalog # M03637. Tested in ELISA, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 3G9H7 / 4H12E9
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q96IZ0


Technical Details

Immunogen	Purified recombinant fragment of PAR4 (aa1-330) expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

Anti-PAR4 PAWR Monoclonal Antibody (M03637) Images



Western Blot analysis using PAR4 Monoclonal Antibody against full-length Trx-PAR4 recombinant protein (1) and HeLa cell lysate (2).



Immunohistochemistry analysis of paraffin-embedded human brain (left) and breast carcinoma (right), showing cytoplasmic and membrane localization with DAB staining using PAR4 Monoclonal Antibody.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PAR4 PAWR Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.