

## Anti-VGLUT1 Rabbit Monoclonal Antibody

Catalog Number: M03664

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



## Overview


|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-VGLUT1 Rabbit Monoclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                             |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Description          | Boster Bio Anti-VGLUT1 Rabbit Monoclonal Antibody catalog # M03664. Tested in WB application. This antibody reacts with Human, Mouse, Rat.                                                                                                                                                                                                                                                                                                         |
| Application          | WB                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Clonality            | Monoclonal 18S50                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Formulation          | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.<br>*This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.                                                                                                                                                                                                                                                                                                         |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Uniprot ID           | Q9P2U7                                                                                                                                                                                                                                                                                                                                                                                                                                             |

## Technical Details

|                     |                                                 |
|---------------------|-------------------------------------------------|
| Immunogen           | A synthesized peptide derived from human VGLUT1 |
| Isotype             | IgG                                             |
| Form                | Liquid                                          |
| Concentration       | 0.5mg/ml                                        |
| Purification        | Affinity-chromatography                         |
| Suggested Dilutions | WB 1:500-1:2000                                 |

## Anti-VGLUT1 Rabbit Monoclonal Antibody (M03664) Images

---



Western blot analysis of VGLUT1 expression in Neuro-2a cell lysate.

### Anti-VGLUT1 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.