

## Anti-HLTF Rabbit Monoclonal Antibody

Catalog Number: M03682

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



## Overview

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-HLTF Rabbit Monoclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                               |
| Reactive Species     | Human, Rat                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Description          | Boster Bio Anti-HLTF Rabbit Monoclonal Antibody catalog # M03682. Tested in WB, IHC, ICC/IF applications. This antibody reacts with Human, Rat.                                                                                                                                                                                                                                                                                                    |
| Application          | IF, IHC, ICC, WB                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Clonality            | Monoclonal 23H05                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Formulation          | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.<br>*This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.                                                                                                                                                                                                                                                                                                         |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Uniprot ID           | Q14527                                                                                                                                                                                                                                                                                                                                                                                                                                             |


## Technical Details

|                     |                                               |
|---------------------|-----------------------------------------------|
| Immunogen           | A synthesized peptide derived from human HLTF |
| Isotype             | IgG                                           |
| Form                | Liquid                                        |
| Concentration       | 0.5mg/ml                                      |
| Purification        | Affinity-chromatography                       |
| Suggested Dilutions | WB 1:1000-5000                                |

|  |                                 |
|--|---------------------------------|
|  | IHC 1:50-200<br>ICC/IF 1:50-200 |
|--|---------------------------------|

## Anti-HLTF Rabbit Monoclonal Antibody (M03682) Images

---



Western blot analysis of HLTF expression in K562 cell lysate.

Anti-HLTF Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.