

Anti-HNRNPU Antibody (C-term)

Catalog Number: M03691-1

About HNRNPU

Component of the CRD-mediated complex that promotes MYC mRNA stabilization. Binds to pre-mRNA. Has high affinity for scaffold-attached region (SAR) DNA. Binds to double- and single- stranded DNA and RNA. Plays a role in the circadian regulation of the core clock component ARNTL/BMAL1 transcription (By similarity).


Overview

Product Name	Anti-HNRNPU Antibody (C-term)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-HNRNPU Antibody (C-term) (Catalog # M03691-1). Tested in WB application(s). This antibody reacts with Human, Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q00839

Technical Details

Immunogen	This HNRNPU antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 766-800 amino acids from the C-terminal region of human HNRNPU.
Predicted Reactive Species	Bovine, Rabbit
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:2000

Anti-HNRNPU Antibody (C-term) (M03691-1) Images



All lanes : Anti-HNRNPU Antibody (C-term) at 1:2000 dilution
Lane 1: Jurkat whole cell lysates
Lane 2: NIH/3T3 whole cell lysates
Lane 3: A431 whole cell lysates
Lysates/proteins at 20 µg per lane. Secondary Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/10000 dilution
Predicted band size : 91 kDa
Blocking/Dilution buffer: 5% NFDm/TBST.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HNRNPU Antibody (C-term)

For Research Use Only. Not for use in diagnostic procedures.