

Anti-GSTA4 Rabbit Monoclonal Antibody

Catalog Number: M03699-1


Overview

Product Name	Anti-GSTA4 Rabbit Monoclonal Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-GSTA4 Rabbit Monoclonal Antibody catalog # M03699-1. Tested in WB applications. This antibody reacts with Human, Rat.
Application	WB
Clonality	Monoclonal 31G31
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O15217


Technical Details

Immunogen	A synthesized peptide derived from human GSTA4
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000

Anti-GSTA4 Rabbit Monoclonal Antibody (M03699-1) Images



Western blot analysis of GSTA4 expression in 293T cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GSTA4 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.