

Anti-ALDH7A1 Rabbit Monoclonal Antibody

Catalog Number: M03712-1


Overview

Product Name	Anti-ALDH7A1 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-ALDH7A1 Rabbit Monoclonal Antibody catalog # M03712-1. Tested in WB, IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal 28A63
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P49419

Technical Details

Immunogen	A synthesized peptide derived from human ALDH7A1
Isotype	IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000 IHC 1:50-1:200

Anti-ALDH7A1 Rabbit Monoclonal Antibody (M03712-1) Images



Western blot analysis of ALDH7A1 expression in SW480 cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ALDH7A1 Rabbit Monoclonal Antibody