

## Anti-Annexin 6 monoclonal ANXA6 Antibody

Catalog Number: M03735


### Overview

Product Name	Anti-Annexin 6 monoclonal ANXA6 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Annexin 6 monoclonal ANXA6 Antibody catalog # M03735. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal AS1A19
Formulation	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P08133

### Technical Details

Immunogen	Recombinant full length Human Annexin 6.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000

## Anti-Annexin 6 monoclonal ANXA6 Antibody (M03735) Images



Western blot detection of ANXA6 in T47D, MB4, CEM, MOLT-4 & Jurkat cell lysates and using ANXA6 antibody (1:1000 diluted) .

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Annexin 6 monoclonal ANXA6 Antibody

For Research Use Only. Not for use in diagnostic procedures.