

Anti-p3 beta MAPK11 Monoclonal Antibody

Catalog Number: M03738


Overview

Product Name	Anti-p3 beta MAPK11 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-p3 beta MAPK11 Monoclonal Antibody catalog # M03738. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 4H6H6
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q15759

Technical Details

Immunogen	Purified recombinant fragment of p38beta (aa251-363) expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-p3 beta MAPK11 Monoclonal Antibody (M03738) Images



Western Blot analysis using p38beta Monoclonal Antibody against truncated p38beta recombinant protein (1) and full-length p38beta (aa1-363)-pcDNA3.1 transfected CHO-K1 cell lysate (2).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p3 beta MAPK11 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.