

Anti-Stabilin 2 (STAB2) Mouse Monoclonal Antibody [Clone ID: OTI2A6]

Catalog Number: M03742-1


Overview

Product Name	Anti-Stabilin 2 (STAB2) Mouse Monoclonal Antibody [Clone ID: OTI2A6]
Reactive Species	Human
Description	Boster Bio STAB2 mouse monoclonal antibody, clone OTI2A6 (formerly 2A6). Catalog# M03742-1. Tested in IHC. This antibody reacts with Human.
Application	IHC
Clonality	Monoclonal OTI2A6
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q8WWQ8

Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 1914-2194 of human STAB2 (NP_060034) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	IHC 1:150

Anti-Stabilin 2 (STAB2) Mouse Monoclonal Antibody [Clone ID: OTI2A6] (M03742-1) Images



Immunohistochemical staining of paraffin-embedded Human endometrium tissue within the normal limits using anti-STAB2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Stabilin 2 (STAB2) Mouse Monoclonal Antibody [Clone ID: OTI2A6]

For Research Use Only. Not for use in diagnostic procedures.