

## Anti-BAP31 Antibody

Catalog Number: M03767

### About BCAP31

Functions as a chaperone protein. Is one of the most abundant endoplasmic reticulum (ER) proteins. Plays a role in the export of secreted proteins in the ER, the recognition of abnormally folded protein and their targeting to the ER associated-degradation (ERAD). Also serves as a cargo receptor for the export of transmembrane proteins. May be involved in CASP8- mediated apoptosis.


### Overview

Product Name	Anti-BAP31 Antibody
Reactive Species	Human
Description	Boster Bio Anti-BAP31 Antibody (Catalog # M03767). Tested in WB, IHC-P application(s). This antibody reacts with Human.
Application	IHC-P, WB
Clonality	Monoclonal 1502CT208.31.16
Formulation	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Mouse
Uniprot ID	P51572


### Technical Details

Immunogen	This BAP31 antibody is generated from a mouse immunized with a recombinant.
Predicted Reactive Species	Bovine, Rabbit
Isotype	IgG1,k
Purification	This antibody is purified through a protein G column, followed by dialysis against PBS.
Suggested Dilutions	WB: 1:2000 IHC-P: 1:25

## Anti-BAP31 Antibody (M03767) Images



All lanes : Anti-BAP31 Antibody at 1:2000 dilution  
Lane 1: A431 whole cell lysates  
Lane 2: NCCIT whole cell lysates  
Lysates/proteins at 20 ug per lane. Secondary Goat Anti-mouse IgG, (H+L), Peroxidase conjugated at 1/10000 dilution  
Predicted band size : 28 kDa  
Blocking/Dilution buffer: 5% NFDm/TBST.



M03767 staining BAP31 in human brain sections by Immunohistochemistry (IHC-P -paraformaldehyde-fixed, paraffin-embedded sections). Tissue was fixed with formaldehyde and blocked with 3% BSA for 0.5 hour at room temperature; antigen retrieval was by heat mediation with a citrate buffer (pH6). Samples were incubated with primary antibody (1/25) for 1 hours at 37°C. A undiluted biotinylated goat polyvalent antibody was used as the secondary antibody.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-BAP31 Antibody

For Research Use Only. Not for use in diagnostic procedures.