

Anti-Aurora C (AURKC) Mouse Monoclonal Antibody [Clone ID: OTI10A7]

Catalog Number: M03809


Overview

Product Name	Anti-Aurora C (AURKC) Mouse Monoclonal Antibody [Clone ID: OTI10A7]
Reactive Species	Human
Description	Boster Bio Anti-AURKC mouse monoclonal antibody, clone OTI10A7 (formerly 10A7). Catalog# M03809. Tested in IF, WB. This antibody reacts with Human.
Application	IF, WB
Clonality	Monoclonal OTI10A7
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9UQB9


Technical Details

Immunogen	Full-length protein expressed in 293T cell transfected with human AURKC expression vector
Isotype	IgG1
Concentration	0.85 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500 IF 1:100

Anti-Aurora C (AURKC) Mouse Monoclonal Antibody [Clone ID: OTI10A7] (M03809) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY AURKC (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-AURKC.



Anti-AURKC mouse monoclonal antibody (M03809) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY AURKC.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Aurora C (AURKC) Mouse Monoclonal Antibody [Clone ID: OTI10A7]

For Research Use Only. Not for use in diagnostic procedures.