

Anti-DEPDC6/DEPTOR Antibody (Monoclonal, 32D55)

Catalog Number: M03811

About DEPTOR

DEP domain-containing mTOR-interacting protein (DEPTOR) also known as DEP domain-containing protein 6 (DEPDC6) is a protein that in humans is encoded by the DEPTOR gene. Enables phosphatidic acid binding activity and protein serine/threonine kinase inhibitor activity. Involved in several processes, including negative regulation of TOR signaling; negative regulation of cell size; and negative regulation of protein kinase activity. Is active in lysosomal membrane.

Overview

Product Name	Anti-DEPDC6/DEPTOR Antibody (Monoclonal, 32D55)
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DEPDC6/DEPTOR Antibody (Monoclonal, 32D55) catalog # M03811. Tested in WB, IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Monoclonal 32D55
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8TB45

Technical Details

Immunogen	Recombinant protein within human DPTOR aa 3-298.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DEPDC6/DEPTOR Antibody (Monoclonal, 32D55)

For Research Use Only. Not for use in diagnostic procedures.