

Anti-COUP TF1 NR2F1 Monoclonal Antibody

Catalog Number: M03833

About NR2F1

Plays a major role in tight junction-specific obliteration of the intercellular space.

Overview


Product Name	Anti-COUP TF1 NR2F1 Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-COUP TF1 NR2F1 Monoclonal Antibody catalog # M03833. Tested in WB, IHC, ICC/IF, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Monoclonal ACHO-14
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P10589

Technical Details

Immunogen	A synthesized peptide derived from human COUP TF1
Isotype	Rabbit IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>WB 1:500-1:1000 IHC 1:50-1:200 ICC/IF 1:50-1:200</p>

FC 1:50

Anti-COUP TF1 NR2F1 Monoclonal Antibody (M03833) Images



Western blot analysis of COUP TF1 expression in HEK293 cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-COUP TF1 NR2F1 Monoclonal Antibody