

Anti-INHBB / Activin B Reference Antibody (Oxford Brookes U. patent anti-Activin Beta-B)

Catalog Number: M03838-1


Overview

Product Name	Anti-INHBB / Activin B Reference Antibody (Oxford Brookes U. patent anti-Activin Beta-B)
Reactive Species	Cynomolgus monkey, Human, Mouse
Description	Boster Bio Anti-INHBB / Activin B Reference Antibody (Oxford Brookes U. patent anti-Activin Beta-B) (Catalog # M03838-1). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human, Mouse. Endotoxin: < 0.825EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P09529


Technical Details

Isotype	IgG1
Form	Liquid

Anti-INHBB / Activin B Reference Antibody (Oxford Brookes U. patent anti-Activin Beta-B) (M03838-1) Images



Anti-INHBB/Activin B Reference Antibody (Oxford Brookes U. patent anti-Activin Beta-B) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-INHBB/Activin B Reference Antibody (Oxford Brookes U. patent anti-Activin Beta-B) is more than 95%

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-INHBB / Activin B Reference Antibody (Oxford Brookes U. patent anti-Activin Beta-B)

For Research Use Only. Not for use in diagnostic procedures.