

Anti-COL5A2 Mouse Monoclonal Antibody [Clone ID: OTI2H3]

Catalog Number: M03869


Overview

Product Name	Anti-COL5A2 Mouse Monoclonal Antibody [Clone ID: OTI2H3]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio COL5A2 mouse monoclonal antibody, clone OTI2H3. Catalog# M03869. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal OTI2H3
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P05997

Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 1200-1499 of human COL5A2 (NP_000384) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/

Anti-COL5A2 Mouse Monoclonal Antibody [Clone ID: OTI2H3] (M03869) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY COL5A2 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-COL5A2 (1:2000).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-COL5A2 Mouse Monoclonal Antibody [Clone ID: OTI2H3]