

## Anti-Arp2 Rabbit Monoclonal Antibody

Catalog Number: M03898


### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-Arp2 Rabbit Monoclonal Antibody   |
| Reactive Species     | Human, Mouse, Rat  |
| Description          | Boster Bio Anti-Arp2 Rabbit Monoclonal Antibody catalog # M03898. Tested in WB, IHC, ICC/IF applications. This antibody reacts with Human, Mouse, Rat. |
| Application          | IF, IHC, ICC, WB   |
| Clonality            | Monoclonal 24A97   |
| Formulation          | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.                                    |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.             |
| Host                 | Rabbit   |
| Uniprot ID           | P61160   |

### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | A synthesized peptide derived from human Arp2  |
| Isotype             | IgG  |
| Form                | Liquid   |
| Concentration       | Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure. |
| Purification        | Affinity-chromatography  |
| Suggested Dilutions | WB 1:500-1:2000<br>IHC 1:50-1:200<br>ICC/IF 1:50-1:200                                       |

## Anti-Arp2 Rabbit Monoclonal Antibody (M03898) Images



Western blot analysis of Arp2 expression in HUVEC cell lysate.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Arp2 Rabbit Monoclonal Antibody