

Anti-TXNRD2 monoclonal Antibody

Catalog Number: M03900


Overview

Product Name	Anti-TXNRD2 monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-TXNRD2 monoclonal Antibody catalog # M03900. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 3B4
Formulation	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9NNW7

Technical Details

Immunogen	Recombinant full length Human TXNRD2.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000

Anti-TXNRD2 monoclonal Antibody (M03900) Images



Western blot detection of TXNRD2 in A431, U251, U2OS, LNCAP&HeLa cell lysates and using TXNRD2 antibody (1:1000 diluted) .

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TXNRD2 monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.