

Anti-Ephrin type-A receptor 7 EphA7 Monoclonal Antibody

Catalog Number: M03943


Overview

Product Name	Anti-Ephrin type-A receptor 7 EphA7 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Ephrin type-A receptor 7 EphA7 Monoclonal Antibody catalog # M03943. Tested in ELISA, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 6C8G7
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q15375


Technical Details

Immunogen	Purified recombinant fragment of EphA7 (aa27-210) expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

Anti-Ephrin type-A receptor 7 EphA7 Monoclonal Antibody (M03943) Images



Western Blot analysis using EphA7 Monoclonal Antibody against truncated GST-EphA7 recombinant protein (1) and truncated EphA7 (aa25-556)-hlgGfc transfected CHOK1 cell lysate (2).



Immunohistochemistry analysis of paraffin-embedded human Breast tissues with AEC staining using EphA7 Monoclonal Antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ephrin type-A receptor 7 EphA7 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.