

## Anti-EGFL7 Mouse Monoclonal Antibody [Clone ID: OTI2B5]

Catalog Number: M04011


### Overview

Product Name	Anti-EGFL7 Mouse Monoclonal Antibody [Clone ID: OTI2B5]
Reactive Species	Human
Description	Boster Bio EGFL7 mouse monoclonal antibody, clone OTI2B5 (formerly 2B5). Catalog# M04011. Tested in WB. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal OTI2B5
Formulation	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9UHF1

### Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 24-273 of human EGFL7 (NP_057299) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:  For protocols, please visit <a href="https://www.bosterbio.com/protocol-and-troubleshooting/">https://www.bosterbio.com/protocol-and-troubleshooting/</a>

## Anti-EGFL7 Mouse Monoclonal Antibody [Clone ID: OTI2B5] (M04011) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY EGFL7 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-EGFL7 (1:2000).

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-EGFL7 Mouse Monoclonal Antibody [Clone ID: OTI2B5]