

## Anti-Parvalbumin Rabbit Monoclonal Antibody

Catalog Number: M04041-3


### Overview

Product Name	Anti-Parvalbumin Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of Parvalbumin protein. Validated for research with WB,IHC,IF,IP,ELISA applications. Anti Parvalbumin is reactive for Human, Mouse, Rat samples
Application	ELISA, IP, IF, IHC, WB
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	P20472


### Technical Details

Isotype	IgG, Kappa
Purification	Protein A
Suggested Dilutions	IHC 1:200-1:1000 WB 1:2000-1:10000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200


## Anti-Parvalbumin Rabbit Monoclonal Antibody (M04041-3) Images




Rat brain was stained with anti-Parvalbumin rabbit antibody



Human brain was stained with anti-Parvalbumin rabbit antibody



Mouse brain was stained with anti-Parvalbumin rabbit antibody



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-Parvalbumin antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: U-251 MG Lane 2: Mouse brain Lane 3: Rat brain Predicted band size: 12kDa Observed band size: 12kDa

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.