

Anti-Rak FRK Monoclonal Antibody

Catalog Number: M04075


Overview

Product Name	Anti-Rak FRK Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Rak FRK Monoclonal Antibody catalog # M04075. Tested in ELISA, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 1A8H2C12
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P42685


Technical Details

Immunogen	Purified recombinant fragment of Rak (aa2-300) expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

Anti-Rak FRK Monoclonal Antibody (M04075) Images



Western blotting validation for Anti-Rak FRK Monoclonal Antibody M04075 Western Blot (WB) analysis using Rak Monoclonal Antibody against HEK293T cells transfected with the pCMV6-ENTRY control (1) and pCMV6-ENTRY Rak cDNA (2). Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration



Immunohistochemistry validation of FRK using Anti-Rak FRK Monoclonal Antibody (M04075). Immunohistochemistry (IHC) analysis of paraffin-embedded Human Breast tissues with AEC staining using Rak Monoclonal Antibody. For more protocol information of IHC

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Rak FRK Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.