

## Anti-FGF18 Antibody (Monoclonal, 32F35)

Catalog Number: M04181

### About Fgf18

Predicted to enable growth factor activity; type 1 fibroblast growth factor receptor binding activity; and type 2 fibroblast growth factor receptor binding activity. Acts upstream of or within several processes, including positive regulation of chondrocyte differentiation; positive regulation of signal transduction; and skeletal system development. Predicted to be located in nucleolus. Predicted to be active in cytoplasm and extracellular space. Is expressed in several structures, including alimentary system; central nervous system; embryo mesenchyme; musculoskeletal system; and sensory organ. Orthologous to human FGF18 (fibroblast growth factor 18).

### Overview

Product Name	Anti-FGF18 Antibody (Monoclonal, 32F35)
Reactive Species	Mouse
Description	Boster Bio Anti-FGF18 Antibody (Monoclonal, 32F35) catalog # M04181. Tested in WB applications. This antibody reacts with Mouse.
Application	WB
Clonality	Monoclonal 32F35
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O89101

### Technical Details

Immunogen	Recombinant protein within mouse FGF18 aa 3-205.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.
Suggested Dilutions	Western blot, 1:500-2000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FGF18 Antibody (Monoclonal, 32F35)

For Research Use Only. Not for use in diagnostic procedures.