

Anti-TRIM22 Mouse Monoclonal Antibody [Clone ID: OTI1D11]

Catalog Number: M04221


Overview

Product Name	Anti-TRIM22 Mouse Monoclonal Antibody [Clone ID: OTI1D11]
Reactive Species	Human
Description	Boster Bio TRIM22 mouse monoclonal antibody, clone OTI1D11 (formerly 1D11). Catalog# M04221. Tested in IF, IHC, WB. This antibody reacts with Human.
Application	IF, IHC, WB
Clonality	Monoclonal OTI1D11
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q8IYM9


Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 61-406 of human TRIM22 (NP_006065) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500 IHC 1:150 IF 1:100


Anti-TRIM22 Mouse Monoclonal Antibody [Clone ID: OTI1D11] (M04221) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TRIM22 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TRIM22.




Immunohistochemical staining of paraffin-embedded Carcinoma of Human thyroid tissue using anti-TRIM22 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)




Immunohistochemical staining of paraffin-embedded Human endometrium tissue within the normal limits using anti-TRIM22 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Immunohistochemical staining of paraffin-embedded Human tonsil within the normal limits using anti-TRIM22 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Immunohistochemical staining of paraffin-embedded Human colon tissue within the normal limits using anti-TRIM22 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Anti-TRIM22 mouse monoclonal antibody (M04221) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY TRIM22.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRIM22 Mouse Monoclonal Antibody [Clone ID: OTI1D11]

For Research Use Only. Not for use in diagnostic procedures.