

Anti-IL-5Ra/ CD125 Reference Antibody (benralizumab)

Catalog Number: M04243-1


Overview

Product Name	Anti-IL-5Ra/ CD125 Reference Antibody (benralizumab)
Reactive Species	Human, Mouse, Rabbit, Rat
Description	Boster Bio Anti-IL-5Ra/ CD125 Reference Antibody (benralizumab) (Catalog # M04243-1). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Rat, Human, Mouse, Rabbit. Endotoxin: < 0.851EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q01344


Technical Details

Isotype	IgG1
Form	Liquid


Anti-IL-5Ra/ CD125 Reference Antibody (benralizumab) (M04243-1) Images



Anti-IL-5Ra/ CD125 Reference Antibody (benralizumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-IL-5Ra/ CD125 Reference Antibody (benralizumab) is more than 98.82%



Immobilized human IL 5Ralpha His at 2 µg/mL can bind Anti-IL-5Ra/ CD125 Reference Antibody (benralizumab)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-5Ra/ CD125 Reference Antibody (benralizumab)

For Research Use Only. Not for use in diagnostic procedures.