

## Anti-AK2 Rabbit Monoclonal Antibody

Catalog Number: M04253-1


### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-AK2 Rabbit Monoclonal Antibody  |
| Reactive Species     | Human, Mouse, Rat  |
| Description          | Boster Bio Anti-AK2 Rabbit Monoclonal Antibody catalog # M04253-1. Tested in WB, IHC, ICC/IF, IP applications. This antibody reacts with Human, Mouse, Rat.  |
| Application          | IP, IF, IHC, ICC, WB   |
| Clonality            | Monoclonal 29A42   |
| Formulation          | Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.<br>*This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.   |
| Host                 | Rabbit   |
| Uniprot ID           | P54819   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | A synthesized peptide derived from human AK2                |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 0.5mg/ml  |
| Purification        | Affinity-chromatography                                     |
| Suggested Dilutions | WB 1:500-2000<br>IHC 1:50-200<br>ICC/IF 1:50-200<br>IP 1:50 |

## Anti-AK2 Rabbit Monoclonal Antibody (M04253-1) Images



Western blot analysis of AK2 expression in HepG2 cell lysate.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AK2 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.